




SERVICE SOLUTIONS

Refurbished Homefill Pump



Many out-of-warranty Invacare Homefill pumps are scrapped due to inadequate repair options. Oxus has engineered a refurbished Homefill pump that achieves as-new OEM performance specifications, and is offered with a **3 year warranty**; doubling the life of your asset, at a fraction of the cost of a new unit.

The fully-tested Oxus refurbished Homefill pump sub-assembly includes:

- 5 cylinder pump assembly
- High pressure inter-stage lines
- Pressure switch

In addition, Oxus offers:

- Detailed written removal and replacement procedure with YouTube training video.
- Key service components for the balance of the Homefill (see reverse side).




*The Oxus Refurbished Homefill Pump
is used in the
Invacare Homefill Oxygen System*

Oxus America
2046 Brown Rd.
Auburn Hills, MI 48326
Phone: 248-475-0925 | Toll Free: 888-475-1568
sales@oxusamerica.com
www.oxusamerica.com




Oxus is an engineering, manufacturing, and service company specializing in medical devices and gas separation technology. We are a FDA registered, ISO 13485 Certified OEM Manufacturer and have been in business for over 20 years.

Refurbished Homefill Pumps and Parts



RS-00519-F

HPP Refurbished Pump - Fixed Frame



RS-00519-S

HPP Refurbished Pump - Slotted Frame

DESCRIPTION

OXUS PART



Pressure Switch

RS-00550



Burst Disk

RS-00551



O2 Tank Fill
Connector

RS-00555




Output Tube

RS-00556



Input Quick
Disconnect

RS-00564



*Our process is not a simple repair, we refurbish the entire pump.
Contact us today to reduce your asset costs*

Oxus America

2046 Brown Rd.

Auburn Hills, MI 48326

Phone: 248-475-0925 | Toll Free: 888-475-1568

sales@oxusamerica.com

www.oxusamerica.com



Oxus is an engineering, manufacturing, and service company specializing in medical devices and gas separation technology. We are a FDA registered, ISO 13485 Certified OEM Manufacturer and have been in business for over 20 years.